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Instructions for use

Product Descriptions
Atrauman dressings are non-medicated, impregnated dressings for atraumatic wound treatment.

Composition
Support fabric made of hydrophobic polyester tulle, impregnated with a non-medicated ointment based on triglycerides (neutral fats). The ointment contains (INCI-analogous names): Caprylic/Capric/Myristic/Stearyl Triglyceride; Bis-Diglyceryl Polyaclyadipate-2.

Properties and mode of action
The soft, thin support material made of polyester allows close contact with the whole surface of the wound which is necessary for the drainage of wound exudate. The water-repellent polyester filaments and the smooth surface structure of the tulle counteract adhesion to the wound. Atrauman allows for atraumatic and painless dressing changes. Atrauman's non-adhering properties are further enhanced by the ointment mass which the product is impregnated. It has no occlusive effects and avoids causing maceration of the underlying skin. Excessive amounts of exudate are drained away into the secondary dressing provided on top of the Atrauman dressing. By keeping the edges of the wound soft and supple, the non-medicated, non-sensitising, paraffin-free ointment base of Atrauman, allows undisturbed epithelialisation. Because it leaves no residue on removal, it conditions the wound for possible surgical procedures and facilitates wound cleansing.

Intended purpose
Atrauman are single-use, sterile, non-medicated ointment dressings suitable for the treatment of superficial acute and chronic wounds of any type. As a non-adherent wound contact layer, it is particularly suitable for preventing the secondary dressing from sticking to the wound bed and for keeping wound and surrounding skin supple. Due to the neutral properties of the non-medicated ointment mass, Atrauman is especially useful in dermatology as well as for patients with sensitive skin and those sensitive to certain medications. Atrauman is suitable for the treatment of human skin by professional users.

Mode of application
Take the dressing out of the peck pack with both the cover papers in place and, if needed, cut it with sterile scissors to the size of the wound. After removal of one of the cover papers, place this side of the dressing onto the wound and remove the second cover. Place a sterile, absorbent dressing pad over the Atrauman dressing to absorb exudate. Unless otherwise prescribed by the physician or another healthcare professional, apply a new Atrauman dressing at each dressing change. The wear time of the dressing depends on the condition of the wound and the level of exudate. Note: If there are therapeutic reasons for not changing the dressing, the ointment may be absorbed when left in place for a longer period of time. In this case possible wound adhesion which may occur can be eliminated by placing a second Atrauman dressing over the first one.

Contraindications
Do not use Atrauman on patients who may be allergic to any of its ingredients.

Product disposal
To minimise the risk of potential infection hazards, or environmental pollution, disposable components of Atrauman should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Special precautions
In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant women, or nursing mothers, the absence of data to the contrary, on these population groups this dressing should be used with caution following a clinician's recommendation.

Notice to the user and patient: Please report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

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Instructions for use
Remove the dressing with both cover papers in place and, if needed, cut it with sterile scissors to the size of the wound. After removal of one of the cover papers, place this side of the dressing onto the wound and remove the second cover. Place a sterile, absorbent dressing pad over the Atrauman dressing to absorb exudate. Unless otherwise prescribed by the physician or another healthcare professional, apply a new Atrauman dressing at each dressing change. The wear time of the dressing depends on the condition of the wound and the level of exudate. Note: If there are therapeutic reasons for not changing the dressing, the ointment may be absorbed when left in place for a longer period of time. In this case possible wound adhesion which may occur can be eliminated by placing a second Atrauman dressing over the first one.

Product Descriptions
Atrauman dressings are non-medicated, impregnated dressings for atraumatic wound treatment.

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